

---

# EEG algorithm SDK for iOS: Development Guide

February 2, 2016

The NeuroSky® product families consist of hardware and software components for simple integration of this biosensor technology into consumer and industrial end-applications. All products are designed and manufactured to meet consumer thresholds for quality, pricing, and feature sets. NeuroSky sets itself apart by providing building block component solutions that offer friendly synergies with related and complementary technological solutions.

**NO WARRANTIES: THE NEUROSKY PRODUCT FAMILIES AND RELATED DOCUMENTATION IS PROVIDED "AS IS" WITHOUT ANY EXPRESS OR IMPLIED WARRANTY OF ANY KIND INCLUDING WARRANTIES OF MERCHANTABILITY, NON-INFRINGEMENT OF INTELLECTUAL PROPERTY, INCLUDING PATENTS, COPYRIGHTS OR OTHERWISE, OR FITNESS FOR ANY PARTICULAR PURPOSE. IN NO EVENT SHALL NEUROSKY OR ITS SUPPLIERS BE LIABLE FOR ANY DAMAGES WHATSOEVER (INCLUDING, WITHOUT LIMITATION, DAMAGES FOR LOSS OF PROFITS, BUSINESS INTERRUPTION, COST OF REPLACEMENT GOODS OR LOSS OF OR DAMAGE TO INFORMATION) ARISING OUT OF THE USE OF OR INABILITY TO USE THE NEUROSKY PRODUCTS OR DOCUMENTATION PROVIDED, EVEN IF NEUROSKY HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES. , SOME OF THE ABOVE LIMITATIONS MAY NOT APPLY TO YOU BECAUSE SOME JURISDICTIONS PROHIBIT THE EXCLUSION OR LIMITATION OF LIABILITY FOR CONSEQUENTIAL OR INCIDENTAL DAMAGES.**

**USAGE OF THE NEUROSKY PRODUCTS IS SUBJECT OF AN END-USER LICENSE AGREEMENT.**

**“Made for iPod,” “Made for iPhone,” and “Made for iPad” mean that an electronic accessory has been designed to connect specifically to iPod, iPhone, or iPad, respectively, and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards. Please note that the use of this accessory with iPod, iPhone, or iPad may affect wireless performance.**

# Contents

<b>About the iOS SDK</b>	<b>4</b>
EEG Algorithm SDK for iOS Contents . . . . .	4
Application development . . . . .	4
Introduction . . . . .	4
EEG Algorithm Sample Project . . . . .	4
<b>API Documentation</b>	<b>6</b>
Data Types . . . . .	6
SDK Delegate Methods . . . . .	7
stateChanged . . . . .	7
attAlgoIndex . . . . .	7
medAlgoIndex . . . . .	7
eyeBlinkDetect . . . . .	8
bpAlgoIndex . . . . .	8
signalQuality . . . . .	8
SDK Utility Methods . . . . .	8
sharedInstance . . . . .	8
getAlgoVersion . . . . .	9
getSdkVersion . . . . .	9
setAlgorithmTypes . . . . .	9
dataStream . . . . .	10
startProcess . . . . .	11
pauseProcess . . . . .	12
stopProcess . . . . .	12
<b>Applications</b>	<b>13</b>
Application of Attention Algorithm . . . . .	13
Application of Meditation Algorithm . . . . .	13
Application of Eye Blink Detection . . . . .	13
Application of EEG Bandpower Computation . . . . .	14
SDK Operations . . . . .	14
Pause and Resume . . . . .	14
Stop and Start . . . . .	14
<b>Frequently Asked Questions</b>	<b>16</b>

# About the iOS SDK

---

This document will guide you through the process of generating algorithm outputs from different NeuroSky Proprietary Mind Algorithms *using* **NeuroSky EEG Algorithm SDK for iOS** with EEG data collected by NeuroSky Biosensor System (e.g. TGAM module or MindWave Mobile Headset).

This development guide is intended for *iOS application developers* who are already familiar with standard iOS development using **Xcode** and Apple's iOS SDK. If you are not already familiar with developing for iOS, please first visit Apple's web site for instruction and tools to develop iOS apps.

**Important:** .

- Requires iOS 6.0 or later

## EEG Algorithm SDK for iOS Contents

- EEG Algorithm SDK for iOS: Development Guide (this document)
- EEG Algorithm SDK framework: AlgoSdk.framework (compatible with both iOS device and simulator)
- Readme: Readme file
- Algo SDK Sample project

## Application development

### Introduction

We recommend developers to use our [COMM SDK for iOS](#) in their application. Comm SDK reduces the complexity of managing EEG Algorithm SDK connections and handles data stream parsing. With the help of our SDKs, the application could be as simple as passing the data received and parsed by the Comm SDK to specific function call(s) at Algo SDK. Specific EEG algorithm index would then be returned accordingly.

**Important:** .

- NeuroSky Comm SDK can only communicate with one paired device at a time.

### EEG Algorithm Sample Project

For collecting EEG data from NeuroSky Biosensor System (e.g. MindWave Mobile headset) with iOS device, please refer to our [COMM SDK for iOS](#) document.

**Algo SDK Sample** is an sample iOS application using Communication (Comm) SDK to connect to NeuroSky Biosensor System (e.g. MindWave Mobile headset) and Algorithm (Algo) SDK for algorithmic computation for specific NeuroSky Mind algorithm, including Attention, Meditation, Appreciation, Mental Effort and Familiarity.

1. Connecting NeuroSky MindWave Mobile headset with the iOS device
2. On the iOS app development machine (e.g. macbook), double click the Algo SDK Sample Xcode project ("**Algo SDK Sample.xcodeproj**") to launch the project with Xcode
3. Select the iOS device (recommended to use iPad) connecting to the working headset as the target device
4. Update the code signing options in the project target settings
5. Select **Product** —> **Run** to build and install the "Algo SDK Sample" app
6. In the sample app:
  - (a) realtime Attention and Meditation indices would be shown on the progress bar
  - (b) eye blink indicator will flash with yellow when blink is detected

**Important: .**

- The sample project requires XCode 6.0 or later.
- The sample project only demonstrates how to iterate with the EEG Algo SDK.
- **Comm SDK** enclosed in the sample project is **version 1.1.7**. Please make sure you are using the latest stable Comm SDK version and make proper changes on sample project if needed.
- It may not be completely compliant with Apple's guidelines for building deploy-able applications.

# API Documentation

The **EEG Algorithm SDK API Reference** in this section contains descriptions of the classes and protocols available in the EEG Algorithm iOS API.

## Data Types

```

/* EEG data signal quality definitions */
typedef NS_ENUM(NSInteger, NskAlgoSignalQuality) {
    NskAlgoSignalQualityGood,          /* Signal quality is in good level */
    NskAlgoSignalQualityMedium,        /* Signal quality is in medium level */
    NskAlgoSignalQualityPoor,          /* Signal quality is in poor level */
    NskAlgoSignalQualityNotDetected    /* Sensor signal is not detected */
};

/* SDK state definitions */
typedef NS_ENUM(NSInteger, NskAlgoState) {
    NskAlgoStateInitiated = 1,         /* Algo SDK initialized */
    NskAlgoStateRunning,               /* Algo SDK is performing analysis (i.e. startProcess()
invoked) */
    NskAlgoStateCollectingBaselineData, /* Algo SDK is collecting baseline data [RESERVED] */
    NskAlgoStateStop,                  /* Algo SDK stops data analysis/baseline collection
[RESERVED] */
    NskAlgoStatePause,                 /* Algo SDK pauses data analysis */
    NskAlgoStateUninitiated,           /* Algo SDK is uninitialized */
    NskAlgoStateAnalysingBulkData      /* Algo SDK is analysing a bulk of EEG data [RESERVED] */
};

/* SDK state change reason definitions */
typedef NS_ENUM(NSInteger, NskAlgoReason) {
    NskAlgoReasonConfigChanged = 1,    /* RESERVED: SDK configuration changed */
    NskAlgoReasonUserProfileChanged,    /* RESERVED: Active user profile has been changed */
    NskAlgoReasonUserTrigger,          /* User triggers */
    NskAlgoReasonBaselineExpired,       /* RESERVED: Baseline expired */
    NskAlgoReasonNoBaseline,            /* RESERVED: No baseline data collected yet */
    NskAlgoReasonSignalQuality          /* Due to signal quality */
};

/* EEG algorithm type definitions */
typedef NS_ENUM(NSInteger, NskAlgoEegType) {
    NskAlgoEegTypeAtt    = 0x008,      /* Attention */
    NskAlgoEegTypeMed    = 0x010,      /* Meditation */
    NskAlgoEegTypeBlink  = 0x080,      /* Eye Blink Detection */
    NskAlgoEegTypeBP     = 0x800       /* EEG Bandpower */
};

/* EEG data type definitions (data from COMM SDK) */
typedef NS_ENUM(NSInteger, NskAlgoDataType) {
    NskAlgoDataTypeEEG,                 /* Raw EEG data */
    NskAlgoDataTypeAtt,                 /* Attention data */

```

```
NskAlgoDataTypeMed,          /* Meditation data */
NskAlgoDataTypePQ,          /* Poor signal quality data */
NskAlgoDataTypeBulkEEG,     /* Bulk EEG data (must be multiple of 512, i.e. Ns of
continuous GOOD EEG data */
};
```

## SDK Delegate Methods

### stateChanged

EEG Algo SDK state change notification delegate method

```
// Required
- (void) stateChanged: (NskAlgoState)state reason: (NskAlgoReason) reason;
```

**Note:** .

- Developer will always need to check with the SDK state and perform proper GUI handling

### attAlgoIndex

Attention Algorithm index notification delegate method

```
// Optional
- (void) attAlgoIndex: (NSNumber*)att_index;
```

**Note:** .

- Attention algorithm has a fixed output interval of 1 second
- Attention algorithm index ranges from 0 to 100 where higher the attention index, higher the attention level

### medAlgoIndex

Meditation Algorithm index notification delegate method

```
// Optional
- (void) medAlgoIndex: (NSNumber*)med_index;
```

**Note:** .

- Meditation algorithm has a fixed output interval of 1 second
- Meditation algorithm index ranges from 0 to 100 where higher the meditation index, higher the meditation level

## eyeBlinkDetect

Eye blink detection notification delegate method

```
// Optional
- (void) eyeBlinkDetect: (NSNumber*) strength;
```

**Note:** .

- **eyeBlinkDetect** will be invoked when eye blink is detected from the provided EEG data

## bpAlgoIndex

EEG Bandpower delegate method

```
// Optional
- (void) bpAlgoIndex: (NSNumber*) delta theta: (NSNumber*) theta alpha: (NSNumber*) alpha
beta: (NSNumber*) beta gamma: (NSNumber*) gamma;
```

**Note:** .

- EEG bandpower values (in dB) have a fixed output interval of 1 second

## signalQuality

EEG data signal quality notification delegate method

```
// Optional
- (void) signalQuality: (NskAlgoSignalQuality) signalQuality;
```

**Note:** .

- Signal Quality was measured and reported at a fixed output interval of 1 second
- SDK state will be changed from **RUNNING** to **PAUSE** when signal quality is poor or sensor off-head is detected. It would return to its previous state (e.g. **RUNNING**) when the signal quality returns to normal

## SDK Utility Methods

### sharedInstance

Developer could always use this method to get the EEG Algorithm instance throughout the app

```
/*
 * Return: EEG Algo instance
 */
+ (id) sharedInstance;
```

**Example**

```
NskAlgoSdk *nskAlgoInstance = [NskAlgoSdk sharedInstance];
```

### getAlgoVersion

Get the specific EEG algorithm version

```
/*
 * Return: Specific EEG algorithm version
 */
- (NSString*) getAlgoVersion: (NskAlgoEegType) algoType;
```

#### Example

```
NSString *version = [NSString stringWithFormat:@"Blink Detection Ver.: %@", [[NskAlgoSdk sharedInstance] getAlgoVersion:NskAlgoEegTypeBlink]];
UIAlertView *alert = [[UIAlertView alloc] initWithTitle:@""
                                                    message:version
                                                    delegate:nil
                                                    cancelButtonTitle:@"OK"
                                                    otherButtonTitles:nil];

[alert show];
```

### getSdkVersion

Get Algo SDK version

```
/*
 * Return: Algo SDK version
 */
- (NSString*) getSdkVersion;
```

#### Example

```
NSString *version = [NSString stringWithFormat:@"SDK Ver.: %@", [[NskAlgoSdk sharedInstance] getSdkVersion]];
UIAlertView *alert = [[UIAlertView alloc] initWithTitle:@""
                                                    message:version
                                                    delegate:nil
                                                    cancelButtonTitle:@"OK"
                                                    otherButtonTitles:nil];

[alert show];
```

### setAlgorithmTypes

Select specific EEG algorithm(s)

```
/*
 * Return: 0 - Algo SDK is initialized successfully; Otherwise, something wrong with SDK initialization (please contact with NeuroSky technical support)
 */
- (NSInteger) setAlgorithmTypes: (NskAlgoEegType) algoTypes;
```

#### Example 1 - Single algorithm

```
// getting the Algo SDK instance
NskAlgoSdk *nskAlgo = [NskAlgoSdk sharedInstance];
// setting self as delegate
[nskAlgo setDelegate: self];
// selecting Attention algorithms
[nskAlgo setAlgorithmTypes: NskAlgoEegTypeAtt];
```

### Example 2 - Multiple algorithms

```
// getting the Algo SDK instance
NskAlgoSdk *nskAlgo = [NskAlgoSdk sharedInstance];
// setting self as delegate
[nskAlgo setDelegate: self];
// selecting Attention and Meditation algorithms
[nskAlgo setAlgorithmTypes: NskAlgoEegTypeAtt|NskAlgoEegTypeMed];
```

#### Note: .

- By invoking **setAlgorithmTypes:** with supported EEG algorithm(s), the SDK will always change back to **INITED**

## dataStream

Feed-in realtime (from COMM SDK) or offline (recorded) EEG data to the EEG Algo SDK

#### Important: .

- EEG Algo SDK handles only the following 4 realtime data output from NeuroSky Biosensor System for now:
  - Poor Signal Quality
  - EEG Raw Data
  - Attention
  - Meditation

```
/*
 * Return: TRUE on success; Otherwise, FALSE
 */
- (BOOL) dataStream: (NskAlgoDataType) type data: (int16_t*) data length: (int16_t) length;
```

### Handling realtime EEG data

```
// In the COMM SDK **onDataReceived** delegate method
-(void) onDataReceived: (NSInteger) datatype data: (int) data obj: (NSObject *) obj
deviceType: (DEVICE_TYPE) deviceType {
    if (deviceType != DEVICE_TYPE_MindWaveMobile) {
        return;
    }
    switch (datatype) {
        case MindDataType_CODE_POOR_SIGNAL:
        {
            int16_t poor_signal[1];
```

```

        poor_signal[0] = (int16_t) data;
        [[ NskAlgoSdk sharedInstance] dataStream: NskAlgoDataTypePQ data: poor_signal length: 1];
    }
    break;

    case MindDataType_CODE_RAW:
    {
        int16_t eeg_data[1];
        eeg_data[0] = (int16_t) data;
        [[ NskAlgoSdk sharedInstance] dataStream: NskAlgoDataTypeEEG data: eeg_data length: 1];
    }
    break;

    case MindDataType_CODE_ATTENTION:
    {
        int16_t attention[1];
        attention[0] = (int16_t) data;
        [[ NskAlgoSdk sharedInstance] dataStream: NskAlgoDataTypeAtt data: attention length: 1];
    }
    break;

    case MindDataType_CODE_MEDITATION:
    {
        int16_t meditation[1];
        meditation[0] = (int16_t) data;
        [[ NskAlgoSdk sharedInstance] dataStream: NskAlgoDataTypeMed data: meditation length: 1];
    }
    break;
}

```

### startProcess

Start analysing feed-in EEG data with selected EEG algorithm(s)

```

/*
 * Return: TRUE on success; Otherwise, FALSE
 */
- (BOOL) startProcess;

```

### Example

```

// getting the Algo SDK instance
NskAlgoSdk *nskAlgo = [NskAlgoSdk sharedInstance];
// setting self as delegate
[nskAlgo setDelegate: self];
// selecting Attention algorithm only
[nskAlgo setAlgorithmTypes: NskAlgoEegTypeAtt];
// start analysing EEG data
[nskAlgo startProcess];

```

#### Note: .

- SDK state will only change to **RUNNING** by invoking **startProcess**

## pauseProcess

Pause analysing feed-in EEG data

```
/*
 * Return: TRUE on success; Otherwise, FALSE
 */
- (BOOL) pauseProcess;
```

### Example

```
// User presses PAUSE button on app to pause the data analysis
- (IBAction)pausePress:(id)sender {
    [[NskAlgoSdk sharedInstance] pauseProcess];
}
```

#### Note: .

- SDK state will change to **PAUSE**
- No **algoIndex** delegate method will not be invoked unless **startProcess** method is invoked again

## stopProcess

Stop analysing feed-in EEG data

```
/*
 * Return: TRUE on success; Otherwise, FALSE
 */
- (BOOL) stopProcess;
```

### Example

```
// User presses STOP button on app to stop the data analysis
- (IBAction)stopPress:(id)sender {
    [[NskAlgoSdk sharedInstance] stopProcess];
}
```

#### Note: .

- SDK state will change to **STOP**
- No **algoIndex** delegate method will be invoked unless **startProcess** method is invoked again

# Applications

---

## Application of Attention Algorithm

- Selecting Attention Algorithm by invoking **setAlgorithmTypes:** method
- Starting EEG data analysis by invoking **startProcess:** method
- Attention index will be returned every 1 second when Algo SDK state is **RUNNING**
- Attention index ranges from **0 to 100**. The higher the index, the higher the attention level

**Note:** .

- Attention has a fixed output interval of 1 second, i.e. one new Attention index every second

## Application of Meditation Algorithm

- Selecting Meditation Algorithm by invoking **setAlgorithmTypes:** method
- Starting EEG data analysis by invoking **startProcess:** method
- Meditation index will be returned every 1 second when Algo SDK state is **RUNNING**
- Meditation index ranges from **0 to 100**. The higher the index, the higher the meditation level

**Note:** .

- Meditation has a fixed output interval of 1 second, i.e. one new Meditation index every second

## Application of Eye Blink Detection

- Selecting Eye Blink Detection by invoking **setAlgorithmTypes:** method
- Starting EEG data analysis by invoking **startProcess:** method
- Eye blink strength will be returned once eye blink is detected when Algo SDK state is **RUNNING**

**Note:** .

- No baseline data collection will be needed

## Application of EEG Bandpower Computation

- Selecting EEG bandpower by invoking **setAlgorithmTypes:** method
- Starting EEG data analysis by invoking **startProcess:** method
- EEG bandpower values (in dB) will be returned every 1 second when Algo SDK state is **RUNNING**

## SDK Operations

### Pause and Resume

- Assuming SDK is in **RUNNING** state
- Pausing EEG algorithm data analysis by invoking **pauseProcess:** method
- Resuming EEG algorithm data analysis by invoking **startProcess:** method




Figure 3.1: Time diagram on Pause/Resume SDK

#### Note: .

- There will be no effect on **pauseProcess:** when previous SDK state is not **RUNNING**

### Stop and Start

- Assuming SDK is in **RUNNING** state
- Stopping EEG algorithm data analysis by invoking **stopProcess:** method
- Restart EEG algorithm data analysis by invoking **startProcess:** method




Figure 3.2: Time diagram on Stop/Start SDK

**Note:** .

- There will be no effect on **stopProcess:** when previous SDK state is not **RUNNING**

# Frequently Asked Questions

---